



Performance Qualification

• *Injector and Flow Precision*

• *Instruments and Fluidics*

<i>Instrument Name</i>	<i>Model</i>	<i>Supplier's Name</i>	<i>Serial Number</i>
Pump	HPG-3400RS	Thermo Scientific	8125299
Autosampler	WPS-3000TRS	Thermo Scientific	8125443
Column Oven	TCC-3000SD	Thermo Scientific	06012832
UV Detector	DAD-3000RS	Thermo Scientific	8043573
Chromeleon Datasystem	V. 6.80 SR16 Build 5387	Thermo Scientific	40640

<i>Accessories</i>	<i>Description</i>
Back Pressure Device	Capillary (L:15 m; ID:0,18 mm)
Sample 4	Caffeine in Water 140 µg/ml
Solvent A	Water (HPLC-Grade)

• *Additional Information*

Customer: Customer's Name
 Operator: Operator's Name
 Operator's Jobtitle

Execution Date: 2020-avr-07
 Next Qualification: 2020-oct

• *Limits, Values and Test Results*

	<i>Limit</i>	<i>Observed Value</i>	<i>Result</i>
Injector Precision - Area	<= 0,500 % RSD	0,233 % RSD	Test passed
Flow Precision - Ret. Time	<= 0,100 % RSD	0,034 % RSD	Test passed
	OR <= 0,0200 min SD	0,0005 min SD	

• *Signatures*


	<i>Signature</i>	<i>Date</i>
Submitter / Operator:	{seq.submitOperator.userName}	{seq.submitTime}
Reviewer:	{seq.reviewOperator.userName}	{seq.reviewTime}
Approver (e-sig. only):	{seq.approveOperator.userName}	{seq.approveTime}

• Data for Injector and Flow Precision Test: Volume 5,0 µl

Sample Name	Ret.Time min Caffeine UV_VIS_1	Area mAU*min Caffeine UV_VIS_1
Injector and flow reproducibility_1	1,5517	116,540
Injector and flow reproducibility_2	1,5533	116,696
Injector and flow reproducibility_3	1,5533	116,385
Injector and flow reproducibility_4	1,5533	116,573
Injector and flow reproducibility_5	1,5533	116,347
Injector and flow reproducibility_6	1,5533	116,787
Injector and flow reproducibility_7	1,5533	116,879
Injector and flow reproducibility_8	1,5533	116,199
Injector and flow reproducibility_9	1,5533	116,235
Injector and flow reproducibility_10	1,5533	116,048
Average:	1,5532	116,469
RSD:	0,034 %	0,233 %
RSD Limit:	<= 0,100 %	<= 0,500 %
SD:	0,0005	
SD Limit:	<= 0,0200	
Result:	ok	ok

- ***Charts for Injector and Flow Precision Test***

• **Chromatogram Overlay**



Injection and Retention Time Precision

Name	Area	Ret.Time	Sample No.	Area
------	------	----------	------------	------

	<i>mAU*min Caffeine UV_VIS_1</i>	<i>min Caffeine UV_VIS_1</i>		<i>Dev. from mean [%]</i>
Injector and flow reproducibility_1	116,540	1,552	1	0,061
Injector and flow reproducibility_2	116,696	1,553	2	0,195
Injector and flow reproducibility_3	116,385	1,553	3	-0,072
Injector and flow reproducibility_4	116,573	1,553	4	0,089
Injector and flow reproducibility_5	116,347	1,553	5	-0,105
Injector and flow reproducibility_6	116,787	1,553	6	0,273
Injector and flow reproducibility_7	116,879	1,553	7	0,352
Injector and flow reproducibility_8	116,199	1,553	8	-0,232
Injector and flow reproducibility_9	116,235	1,553	9	-0,201
Injector and flow reproducibility_10	116,048	1,553	10	-0,362
Average:	116,469	1,553		
RSD:	0,233 %	0,034 %		
SD:		0,001 min		

Calculation of RSD and Average

	Area:	Retention Time:
Sum of x:	1164,688041	15,53166667
sum of sq(x):	135650,4876	24,12326944
Sq of sum (x):	1356498,234	241,2326694
STD:	0,27166761	0,000527046
Average:	116,4688041	1,553166667
RSD:	0,002332535	0,000339337
Number of samples:	10	10

Determination of Pump Unit for Dionex DGPs

Sequence name: PQ_INJECTOR_FLOW_REPRO
 Right end of the sequence name : EPRO
 Pump's model number: HPG-3400RS
 Pump's model variant: HPG
 Pump unit:

Evaluation of Accessories (Following autosampler uses standard 3 or 5 instead of st

ModelNo	Loop Volume	Evaluation*)	Reference value for report
WPS-3000	40	0	S4
ACC-3000	200	0	
ACC-3000T	200	0	
WPS-3000	344	0	
WPS-3000	250	0	
Gina 50		0	
Gina 160		0	
Accela Sampler		0	
Accela Open Sampler		0	
OAS-3000		0	
VH-A10-A		0	
VR-A12-A		0	
VR-A13-A		0	
VR-A14-A		0	

VR-A15-A	0
VF-A10-A	0

*) ModelNos without loop volume setting uses different evaluation formulas compared to ModelNos with loop volume setting.









No.	Ret. Time
-----	-----------

	<i>Deviation from the mean value [%]</i>
1	-0,097
2	0,011
3	0,011
4	0,011
5	0,011
6	0,011
7	0,011
8	0,011
9	0,011
10	0,011

Standard 4):